

## **Stress measurement using psychometrics & AR/VR Simulators**

Stress and anxiety levels **can be measured indirectly** using various physiological sensors that track bodily responses correlated with stress. While these measures aren't 100% definitive for diagnosing anxiety (which is complex and often psychological), **they can offer reliable indicators**. To enhance the accuracy and reliability of stress measurement a hybrid system can be so developed that it has the following sub-systems :

1. Psychological set-up to expose individual to different emotional situations
2. AR/VR set-up to display the situations are videos/images
3. Physiological measurement device to measure various changes in physiology of individual
4. Automation set-up to collect all information and project it in desired format to indicate stress levels and vulnerability of an individual for self-harm and harm to others.

Following is a tentative breakdown of the **methods, equipment**, and some **AR/VR examples**:

### 1. Physiological Markers and Sensors for Stress/Anxiety

<u>Biomarker</u>	<u>Measurement</u>	<u>Associated Equipment</u>
<b>Heart Rate (HR)</b>	Elevated HR is a common stress response	ECG, PPG sensors, smartwatches, chest straps
<b>Heart Rate Variability (HRV)</b>	Lower HRV often correlates with high stress	ECG, chest straps, smart rings (e.g., Oura Ring)
<b>Electrodermal Activity (EDA/GSR)</b>	Skin conductance increases with stress	GSR sensors (e.g., Empatica E4, Shimmer)
<b>Blood Pressure (BP)</b>	Often elevated during stress	Digital BP cuffs, wearable BP monitors
<b>Pulse Oximetry (SpO2)</b>	Mild drop in oxygen saturation possible under extreme stress	Pulse oximeters
<b>Electroencephalography (EEG)</b>	Certain brainwave patterns (e.g., high beta waves) can indicate anxiety	EEG headsets (e.g., Muse, Emotiv)
<b>Respiratory Rate</b>	Rapid breathing is a stress response	Respiratory belts, smart shirts (e.g., Hexoskin)
<b>Skin Temperature</b>	Drops in peripheral skin temperature may occur	Infrared sensors, thermistors

## **Measurement Methods**

### **1. HRV Analysis**

- **Time-domain:** RMSSD (root mean square of successive differences)
- **Frequency-domain:** LF/HF ratio
- Devices: Polar H10 chest strap, Garmin smartwatches, Oura Ring

### **2. EEG Analysis**

- Brainwave patterns (e.g., increase in beta, decrease in alpha or theta)
- Devices: Muse 2, Emotiv Insight, NeuroSky MindWave

### **3. EDA Measurement**

- Skin conductance response (SCR) spikes indicate arousal
- Devices: Empatica E4, Affectiva Q Sensor

### **4. Multimodal Approaches**

- Combine HRV + GSR + EEG + respiration for higher accuracy
- Often used in research and clinical settings

## **2. AR/VR in Stress/Anxiety Monitoring and Management**

### **Use of AR/VR for Biofeedback and Anxiety Reduction**

AR/VR environments provide immersive experiences that can be combined with biosensors for real-time monitoring and **adaptive feedback**.

#### **Examples of Equipment & Applications:**

<b><u>System</u></b>	<b><u>Tech Used</u></b>	<b><u>Functionality</u></b>
<b>Brain activity</b>	EEG + VR	Monitors brain activity and provides immersive relaxation/training
<b>General Physiology</b>	VR + physiological sensors	Used in hospitals for anxiety/pain reduction during procedures
<b>Heart</b>	VR + heart rate sensor	Clinical anxiety treatment using biofeedback and exposure therapy
<b>Heart rate, brain</b>	VR + biofeedback (HRV, EEG)	VR games that adapt based on user's stress response
<b>Breathing</b>	VR + breathing sensor	Underwater VR game that uses breath control to navigate – encourages relaxation

**Conceptual system architecture / flowchart for a Stress/Anxiety Monitoring and Management System using Wearable Sensors + VR/AR + Biofeedback Loop:**

1. **User wears sensors and uses a VR/AR device.**
2. **User is exposed to 6 to 7 different situations (created by professional psychologists).**
3. **User observes the situations and his body/mind react to it.**
4. **Wearable Sensor Layer** collects biometric data such as:
  - Heart Rate (via PPG or ECG)
  - Heart Rate Variability (HRV)
  - Galvanic Skin Response (GSR) / Electrodermal Activity (EDA)
  - EEG (brain wave activity)
  - Optional: Blood Oxygen (SpO<sub>2</sub>) and Blood Pressure (BP)
5. **Signal Processing Module** cleans and prepares the data:
  - Filters out noise
  - Extracts relevant features
  - Syncs signals in real-time
6. **Stress/Anxiety Classifier** analyzes the data using either:
  - A machine learning model (designed by expert psychologists and physicians) or
  - A rule-based system {designed by expert psychologists and physicians} (e.g., high heart rate + low HRV = likely stress)

- Also uses EEG, GSR, BP data for better accuracy

7. **Feedback Controller & Recommender** takes action based on detected stress levels:

- Adjusts the VR/AR experience to calm the user
- Recommends breathing exercises
- Triggers guided relaxation content

8. **Immersive VR/AR App** responds dynamically:

- Changes the virtual environment
- Displays calming, stress-relieving content (e.g., nature scenes, meditative visuals)

### **Rough financials (ballpark breakdown)**

**1. Design & development** : Design of whole system including psychological situations, measurement parameters, selection of measuring equipment, signal processing, analytics model, feedback and control, AR/VR system design : **₹ 35 lacs** (including psychologist's & physician's consultancy)

**2. Wearable Physiological Sensors(clinical grade)** : (Research-grade sensors: HR, HRV, EDA, PPG, skin temp, SPO2, BP etc ) : **₹ 30 lacs** (including import duties)

**3. . VR Headsets & PC / Standalone Unit & PC-VR Setup** : **₹ 10 lacs**

**3. VR Software & Biofeedback Integration** :

**Psious VR Platform (clinical-grade exposure therapy with sensors)** : **₹ 5 lacs/year**

**4. Sensors + Software Integration & Support** : **₹ 5 lacs**

**5. Miscellaneous including training** : **₹ 5 lacs**

**Total : 90 lacs**

Note : We can tailor the system by choosing components aligned with specific needs of the Army.